


Product Information	
Catalog No.	FC10118
Product name	Anti-Human EPHA3 Antibody (SAA0126)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0126
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HEK, Eph-like tyrosine kinase 1, EK4, EPHA3, Tyrosine-protein kinase receptor ETK1, Tyrosine-protein kinase TYRO4, Ephrin type-A receptor 3, hEK4, ETK, TYRO4, ETK1, Human embryo kinase, EPH-like kinase 4
Accession	P29320
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human EPHA3 Antibody (SAA0126)

SDS-PAGE