

Product Information	
Catalog No.	FC10137
Product name	Anti-Human CD122/IL2RB Antibody (MikB1)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	MikB1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.38 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL-2R subunit beta, Interleukin-15 receptor subunit beta, IL2RB, IL-2 receptor subunit beta, CD122, High affinity IL-2 receptor subunit beta, p75, Interleukin-2 receptor subunit beta, IL-2RB, p70-75, IL15RB
Accession	P14784
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.