

Product Information	
Catalog No.	FC10170
Product name	Anti-Human CD152/CTLA4 Antibody (SAA0062)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0062
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.38 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD152, Cytotoxic T-lymphocyte protein 4, CTLA-4, CTLA4, Cytotoxic T-lymphocyte-associated antigen 4
Accession	P16410
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.