


Product Information	
Catalog No.	FC10172
Product name	Anti-Human CD154/CD40LG/TNFSF5 Antibody (SAA0549)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0549
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	sCD40L, Tumor necrosis factor ligand superfamily member 5, T-cell antigen Gp39, CD40L, TNFSF5, CD154, CD40-L, TNF-related activation protein, CD40 ligand, CD40LG, TRAP
Accession	P29965
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD154/CD40LG/TNFSF5 Antibody (SAA0549)

SDS-PAGE