

Product Information	
Catalog No.	FC10192
Product name	Anti-Human CD178/FASLG/TNFSF6 Antibody (NOK1/NOK-1)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	NOK1/NOK-1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Receptor-binding FasL ectodomain, APL, APT1LG1, FasL ICD, Fas antigen ligand, CD178, sFasL, CD95L, APTL, TNFSF6, Apoptosis antigen ligand, SPA, FASLG, FasL, SPPL2A-processed FasL form, Tumor necrosis factor ligand superfamily member 6, Soluble Fas ligand, Fas ligand, CD95 ligand, CD95-L, FASL
Accession	P48023
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.