

Product Information	
Catalog No.	FC10230
Product name	Anti-Human CD212/IL12RB1 Antibody (148C09)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH-mFc
Clone ID	148C09
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL-12 receptor beta component, CD212, IL-12R-beta-1, IL-12RB1, IL12RB1, Interleukin-12 receptor subunit beta-1, IL-12 receptor subunit beta-1, IL12R, IL-12R subunit beta-1, IL12RB
Accession	P42701
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.