


Product Information	
Catalog No.	FC10231
Product name	Anti-Human CD213a1/IL13RA1 Antibody (SAA0070)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0070
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL13R, IL-13R subunit alpha-1, IL13RA1, IL13RA, IL-13RA1, CD213a1, Interleukin-13 receptor subunit alpha-1, IL-13 receptor subunit alpha-1, CT19, IL-13R-alpha-1, Cancer/testis antigen 19
Accession	P78552
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD213a1/IL13RA1 Antibody (SAA0070)

SDS-PAGE