


Product Information	
Catalog No.	FC10257
Product name	Anti-Human CD243/P-Glycoprotein Antibody (C219)
Applications	FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	C219
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Phospholipid transporter ABCB1, MDR1, CD243, ATP-binding cassette sub-family B member 1, ABCB1, Multidrug resistance protein 1, P-glycoprotein 1, ATP-dependent translocase ABCB1, PGY1, P-Glycoprotein, P Glycoprotein
Accession	P08183
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD243/P-Glycoprotein Antibody (C219)

SDS-PAGE