


Product Information	
Catalog No.	FC10285
Product name	Anti-Human CD274/PD-L1/B7-H1 Antibody (SAA0088)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA0088
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.02 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	B7-H1, Programmed cell death 1 ligand 1, PDCD1 ligand 1, PDCD1L1, B7 homolog 1, PDCD1LG1, PDL1, hPD-L1, Programmed death ligand 1, B7H1, PD-L1, CD274
Accession	Q9NZQ7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD274/PD-L1/B7-H1 Antibody (SAA0088)

SDS-PAGE