

Product Information	
Catalog No.	FC10300
Product name	Anti-Human CD29/ITGB1 Antibody (K20)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	K20
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	MSK12, VLA-4 subunit beta, CD29, Integrin beta-1, ITGB1, FN1R, MD2, GPIIA, Glycoprotein IIA, Fibronectin receptor subunit beta
Accession	P05556
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.