

Product Information	
Catalog No.	FC10342
Product name	Anti-Human CD350/FZD10 Antibody (ONC-005)
Applications	ELISA, FCM, IHC, IP, RIA
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	ONC-005
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	FZD10, Fz-10, hFz10, Frizzled-10, CD350, FzE7
Accession	Q9ULW2
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.