

Product Information	
Catalog No.	FC10353
Product name	Anti-Human CD367/CLEC4A Antibody (9E8#)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	9E8#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	C-type lectin DDB27, CLECSF6, CLEC4A, LLIR, CD367, Lectin-like immunoreceptor, DCIR, C-type lectin domain family 4 member A, Dendritic cell immunoreceptor, C-type lectin superfamily member 6
Accession	Q9UMR7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.