


Product Information


Catalog No.	FC10357
Product name	Anti-Human CD37/TSPAN26 Antibody (SAA0012)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0012
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD37, Tspan-26, TSPAN26, Leukocyte antigen CD37, Tetraspanin-26
Accession	P11049
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



Flow Cytometry

Flow-cytometry using anti-human CD37 antibody. Human RPMI-8226 cell line were stained with an irrelevant antibody (Green Histogram) or an anti-human CD37 antibody monoclonal antibody (Catalog # HY564207 , Red Histogram) at a concentration of 20 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a FITC conjugated goat anti-human antibody (Catalog # HF690414) and cells analysed on a NovoCyte Flow Cytometer.n



Western blot

Western blot analysis was performed using anti-IL23A monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant human IL23A (Catalog No: HV466011)
Lane 2: negative control