


### Product Information

|                       |                                                                                                                                                                                                                              |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC10420                                                                                                                                                                                                                      |
| Product name          | Anti-Human CD69 Antibody (SAA0031)                                                                                                                                                                                           |
| Applications          | FCM                                                                                                                                                                                                                          |
| Species reactivity    | Human                                                                                                                                                                                                                        |
| Host species          | Mouse                                                                                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                                                                                   |
| Isotype               | IgG1                                                                                                                                                                                                                         |
| Clone ID              | SAA0031                                                                                                                                                                                                                      |
| Form                  | Liquid                                                                                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                      |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                                                                                              |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                                                                                          |
| Target                | AIM, CD69, C-type lectin domain family 2 member C, Early T-cell activation antigen p60, BL-AC/P26, GP32/28, MLR-3, Early activation antigen CD69, Leukocyte surface antigen Leu-23, CLEC2C, EA1, Activation inducer molecule |
| Accession             | Q07108                                                                                                                                                                                                                       |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                                                                                            |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.                                                                 |
| Note                  | For research use only.                                                                                                                                                                                                       |

### Image Data



SEC-HPLC

SEC-HPLC detection for Anti-Human CD69 Antibody (ABS2045).