


Product Information	
Catalog No.	FC10430
Product name	Anti-Human CD74 Antibody (LL1)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	LL1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HLA-DR antigens-associated invariant chain, HLA class II histocompatibility antigen gamma chain, CD74, CLIP, Ia antigen-associated invariant chain, DHLAP, Ii
Accession	P04233
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SEC-HPLC detection for Anti-Human CD74 Antibody (LL1).

Bioactivity