

Product Information	
Catalog No.	FC10441
Product name	Anti-Human CD85k/LILRB4 Antibody (128-3)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	128-3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD85 antigen-like family member K, LIR-5, Leukocyte immunoglobulin-like receptor 5, ILT3, ILT-3, Immunoglobulin-like transcript 3, LILRB4, CD85k, Leukocyte immunoglobulin-like receptor subfamily B member 4, Monocyte inhibitory receptor HM18, LIR5
Accession	Q8NHJ6
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.