

Product Information	
Catalog No.	FC10451
Product name	Anti-Human CD95/FAS Antibody (SAA0040)
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0040
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	APT1, TNFRSF6, FAS1, FASLG receptor, CD95, Apo-1 antigen, Tumor necrosis factor receptor superfamily member 6, FAS, Apoptosis-mediating surface antigen FAS
Accession	P25445
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.