


Product Information	
Catalog No.	FC10460
Product name	Anti-Mouse CD25/IL2RA Antibody (PC61/PC61.5.3)
Applications	FCM, FuncS
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	PC61/PC61.5.3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.01 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, p55, CD25, Il2ra, Il2r
Accession	P01590
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse CD25/IL2RA Antibody (PC61/PC61.5.3)

SDS-PAGE