

Product Information	
Catalog No.	FC10593
Product name	Anti-Mouse CD71/TFRC Antibody (OX-26)
Applications	FCM, IF, IHC (P), IP, WB
Species reactivity	Human, Mouse, Rat
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	OX-26
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.8 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Transferrin receptor protein 1, TR, TfR, TfR1, Trfr, CD71, Tfrc
Accession	Q62351
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.