


Product Information	
Catalog No.	FC10594
Product name	Anti-Mouse CD71/TFRC Antibody (8D3)
Applications	FCM, FuncS, ICC, IF, IHC, IHC-Fr
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	8D3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Transferrin receptor protein 1, TfR1, T9, Trfr, p90, CD71, TFRC, TfR, sTfR, TR
Accession	Q62351
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse CD71/TFRC Antibody (8D3)

SDS-PAGE