


Product Information	
Catalog No.	FC10618
Product name	Anti-Mouse CD207 Antibody (2G3k)
Applications	ELISA, FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	2G3k
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	C-type lectin domain family 4 member K, Langerin, CD207, Clec4k
Accession	Q8VBX4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse CD207 Antibody (2G3k)

SDS-PAGE