


Product Information

Catalog No.	FC10622
Product name	Anti-Mouse CD230 Antibody (POM1)
Applications	ELISA, FCM
Species reactivity	Human, Mouse, Bovine, Hamster
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	POM1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Major prion protein, PrP, PrP27-30, PrP33-35C, CD230, Prnp, Prn-p
Accession	P04925
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-CD230/PRNP monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant human CD230/PRNP (Catalog No: HY392011)
Lane 2: negative control