

Product Information	
Catalog No.	FC10655
Product name	Anti-Mouse ROR1 Antibody (R11)
Applications	ELISA, FCM
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	R11
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.37 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Neurotrophic tyrosine kinase, receptor-related 1, Inactive tyrosine-protein kinase transmembrane receptor ROR1, NTRKR1, ROR1
Accession	Q9Z139
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.