


Product Information	
Catalog No.	FC10661
Product name	Anti-Mouse Ly6g/Ly-6G.1 Antibody (1A8)
Applications	in vivo neutrophil depletion, in vivo MDSC depletion, IHC, IF, FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	1A8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Lymphocyte antigen 6G, Ly-6G, Ly-6G.1, Ly6g, GR-1, GR1
Accession	P35461
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse Ly6g/Ly-6G.1 Antibody (1A8)

SDS-PAGE