

Product Information	
Catalog No.	FC10662
Product name	Anti-Mouse TIMD4/TIM-4 Antibody (RMT4-53)
Applications	FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	RMT4-53
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TIMD4, T-cell immunoglobulin mucin receptor 4, TIM-4, TIMD-4, T-cell membrane protein 4, TIM4, T-cell immunoglobulin and mucin domain-containing protein 4
Accession	Q6U7R4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.