


Product Information	
Catalog No.	FC10678
Product name	Anti-Mouse CD96 Antibody (SAA0311)
Applications	FCM, Neutralization
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0311
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	T-cell surface protein tactile, Cell surface antigen CD96, T cell-activated increased late expression protein, CD96
Accession	Q3U0X8
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse CD96 Antibody (SAA0311)

SDS-PAGE