


Product Information	
Catalog No.	FC10700
Product name	Anti-Mouse NKG2D Antibody (CX5)
Applications	Blocking, ELISA, FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	CX5
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	NKG2-D type II integral membrane protein, Killer cell lectin-like receptor subfamily K member 1, NK cell receptor D, NKG2-D-activating NK receptor, CD314, Klrk1, Nkg2d
Accession	O54709
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mouse NKG2D Antibody (CX5)

SDS-PAGE