


Product Information	
Catalog No.	FC10714
Product name	Anti-Human IL5 Antibody (SAA0374)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA0374
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL5, B-cell differentiation factor I, IL-5, T-cell replacing factor, TRF, Eosinophil differentiation factor, Interleukin-5
Accession	P05113
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human IL5 Antibody (SAA0374)

SDS-PAGE