


Product Information


Catalog No.	FC10722
Product name	Anti-Human IL13 Antibody (SAA0382)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0382
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL13, NC30, Interleukin-13, IL-13
Accession	P35225
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



Flow Cytometry

Flow-cytometry using anti-human IL13 antibody. IL13 Transfected CHO cells were stained with an irrelevant antibody (Blue Histogram) or an anti-human IL13 antibody monoclonal antibody (Catalog # HW632107, Green Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a FITC conjugated goat anti-mouse antibody (Catalog # MF690414) and cells analysed on a NovoCyte Flow Cytometer.



Western blot analysis was performed using anti-IL13 monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant human IL13 (Catalog No: HW632012)
Lane 2: negative control

Western blot