


Product Information	
Catalog No.	FC10744
Product name	Anti-Human CD85j/LILRB1/ILT2 Antibody (SAA0402)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0402
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	LIR1, Immunoglobulin-like transcript 2, MIR7, CD85 antigen-like family member J, Monocyte/macrophage immunoglobulin-like receptor 7, ILT-2, Leukocyte immunoglobulin-like receptor subfamily B member 1, LIR-1, CD85j, ILT2, LILRB1, MIR-7, Leukocyte immunoglobulin-like receptor 1
Accession	Q8NHL6
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD85j/LILRB1/ILT2 Antibody (SAA0402)

SDS-PAGE