


| Product Information   |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC10750                                                                                                                                                      |
| Product name          | Anti-Mouse CXCL10/IP-10 Antibody (SAA0412)                                                                                                                   |
| Applications          | ELISA, FCM                                                                                                                                                   |
| Species reactivity    | Mouse                                                                                                                                                        |
| Host species          | Rat                                                                                                                                                          |
| Clonality             | Monoclonal                                                                                                                                                   |
| Isotype               | IgG1, kappa                                                                                                                                                  |
| Clone ID              | SAA0412                                                                                                                                                      |
| Form                  | Liquid                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Concentration         | 0.75 mg/ml                                                                                                                                                   |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                              |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                | Gamma-IP10, C-X-C motif chemokine 10, INP10, 10 kDa interferon gamma-induced protein, Small-inducible cytokine B10, SCYB10, IP-10, CXCL10                    |
| Accession             | P17515                                                                                                                                                       |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                            |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Note                  | For research use only.                                                                                                                                       |

### Image Data



SDS-PAGE for Anti-Mouse CXCL10/IP-10 Antibody (SAA0412)

SDS-PAGE