


| Product Information   |                                                                                                                                                                                                                 |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC10787                                                                                                                                                                                                         |
| Product name          | Anti-Human PDGFB Antibody (SAA0449)                                                                                                                                                                             |
| Applications          | ELISA, FCM                                                                                                                                                                                                      |
| Species reactivity    | Human                                                                                                                                                                                                           |
| Host species          | Mouse                                                                                                                                                                                                           |
| Clonality             | Monoclonal                                                                                                                                                                                                      |
| Isotype               | IgG2a, kappa                                                                                                                                                                                                    |
| Clone ID              | SAA0449                                                                                                                                                                                                         |
| Form                  | Liquid                                                                                                                                                                                                          |
| Storage buffer        | 0.01M PBS, pH 7.4, 0.09% Sodium azide.                                                                                                                                                                          |
| Concentration         | 1 mg/ml                                                                                                                                                                                                         |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                                                                                 |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                                                                             |
| Target                | PDGF-2, Proto-oncogene c-Sis, PDGF2, SIS, Platelet-derived growth factor B chain, PDGFB, Becaplermin, Platelet-derived growth factor beta polypeptide, Platelet-derived growth factor subunit B, PDGF subunit B |
| Accession             | P01127                                                                                                                                                                                                          |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                                                                               |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.                                                    |
| Note                  | For research use only.                                                                                                                                                                                          |

### Image Data



SDS-PAGE for Anti-Human PDGFB Antibody (SAA0449)

SDS-PAGE