


Product Information	
Catalog No.	FC10850
Product name	Anti-Human SCARB2/CD36L2 Antibody (SAA0777)
Applications	Blocking, ELISA, FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA0777
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CD36, 85 kDa lysosomal membrane sialoglycoprotein, CD36 antigen-like 2, CD36L2, SCARB2, Lysosome membrane protein 2, LIMP II, Scavenger receptor class B member 2, Lysosome membrane protein II, LIMPII, LIMP2, LGP85
Accession	Q14108
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human SCARB2/CD36L2 Antibody (SAA0777)

SDS-PAGE