

Product Information	
Catalog No.	FC10860
Product name	Anti-Human SLITRK6 Antibody (SAA0827)
Applications	ELISA, FCM, IF
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG2, kappa
Clone ID	SAA0827
Form	Liquid
Storage buffer	30 mM Histidine, pH5.8, 10% sucrose, 0.02% Tween80
Concentration	2.14 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	SLITRK6, SLIT and NTRK-like protein 6
Accession	Q9H5Y7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.