


| Product Information   |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | FC10863                                                                                                                                                      |
| Product name          | Anti-Human TNFSF12/APO3L Antibody (SAA0806)                                                                                                                  |
| Applications          | ELISA, FCM                                                                                                                                                   |
| Species reactivity    | Human                                                                                                                                                        |
| Host species          | Mouse                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                   |
| Isotype               | IgG1, kappa                                                                                                                                                  |
| Clone ID              | SAA0806                                                                                                                                                      |
| Form                  | Liquid                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4, 0.09% Sodium azide.                                                                                                                       |
| Concentration         | 1 mg/ml                                                                                                                                                      |
| Purity                | >95% as determined by SDS-PAGE.                                                                                                                              |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                | TWEAK, DR3LG, TNFSF12, Tumor necrosis factor ligand superfamily member 12, APO3L, TNF-related weak inducer of apoptosis, APO3 ligand                         |
| Accession             | O43508                                                                                                                                                       |
| Endotoxin level       | Please contact with the lab for this information.                                                                                                            |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Note                  | For research use only.                                                                                                                                       |

## Image Data

Anti-Human TNFSF12/APO3L Antibody (SAA0806)  
binds with TNFSF12/APO3L



Detects Recombinant Human TNFSF12/APO3L/TWEAK Protein, C-His (Cat No.: HT065011 ) in indirect ELISA.