

Product Information	
Catalog No.	FC10880
Product name	Anti-Human ADRB1 Antibody (6B9)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	VHH-hFc
Clone ID	6B9
Form	Liquid
Storage buffer	9mM Citric acid monohydrate, 41mM Na ₂ HPO ₄ , 1% Mannitol, 5% Trehalose, 0.02% Tween 20, pH 6.5.
Concentration	2.22 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Beta-1 adrenoreceptor, ADRB1, Beta-1 adrenergic receptor, ADRB1R, B1AR, Beta-1 adrenoreceptor
Accession	P08588
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.