


Product Information	
Catalog No.	FC10930
Product name	Anti-Human CD335/NCR1/NKp46 Antibody (8E5#)
Applications	FCM, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	8E5#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	LY94, Natural cytotoxicity triggering receptor 1, Natural killer cell p46-related protein, NK cell-activating receptor, NK-p46, CD335, NKp46, NCR1, hNKp46, Lymphocyte antigen 94 homolog
Accession	O76036
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.


Image Data



SDS-PAGE for Anti-Human CD335/NCR1/NKp46 Antibody (8E5#)


SDS-PAGE

Anti-Human CD335/NCR1/NKp46 Antibody (8E5#)
binds with CD335



Bioactivity

Detects recombinant human CD335/NCR1/NKp46 (Catalog No: HF858021) in indirect ELISA



Western blot

Western blot analysis was performed using anti-CD335/NCR1/NKp46 monoclonal antibody at 1ug/mL on various samples.

Lane 1: recombinant human CD335/NCR1/NKp46 (Catalog No: HF858021)

Lane 2: negative control

Lane 3: negative control