


Product Information	
Catalog No.	FC10932
Product name	Anti-Human HS6ST2 Antibody (SAA1430)
Applications	ELISA, FCM, WB
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1430
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HS6ST-2, Heparan-sulfate 6-O-sulfotransferase 2, HS6ST2
Accession	Q96MM7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human HS6ST2 Antibody (SAA1430)

SDS-PAGE