


Product Information	
Catalog No.	FC10934
Product name	Anti-Human CD15/SSEA-1/Lewis x Antibody (SAA1433)
Applications	ELISA, FCM, IHC
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1433
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.32 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Alpha-(1, 3)-fucosyltransferase 4, ELAM-1 ligand fucosyltransferase, Fucosyltransferase 4, Fucosyltransferase IV, Galactoside 3-L-fucosyltransferase, Fuc-TIV, FCT3A, FUT4, ELFT, FucT-IV
Accession	P22083
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD15/SSEA-1/Lewis x Antibody (SAA1433)

SDS-PAGE