


Product Information	
Catalog No.	FC10959
Product name	Anti-Human CD127/IL7R Antibody (4A10)
Applications	ELISA, FCM, Neutralization
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	4A10
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL7R, Interleukin-7 receptor subunit alpha, IL-7R-alpha, IL-7 receptor subunit alpha, CD127, CDw127, IL-7R subunit alpha, IL-7RA
Accession	P16871
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD127/IL7R Antibody (4A10)

SDS-PAGE