

Product Information	
Catalog No.	FC10981
Product name	Anti-Human SLC39A6 Antibody (SAA1475)
Applications	ELISA, FCM, Neutralization
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA1475
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	SLC39A6, LIV1, Zrt- and Irt-like protein 6, Zinc transporter ZIP6, ZIP-6, Estrogen-regulated protein LIV-1, Solute carrier family 39 member 6, ZIP6
Accession	Q13433
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.