



Product Information	
Catalog No.	FC10986
Product name	Anti-Human MSLN/Mesothelin Antibody (SS1)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SS1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1.21 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CAK1 antigen, Mesothelin, MSLN, Pre-pro-megakaryocyte-potentiating factor, MPF
Accession	Q13421
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human MSLN/Mesothelin Antibody (SS1)

SDS-PAGE




Flow cytometry

Flow-cytometry using anti-human MSLN antibody. Hela cells were stained with an irrelevant antibody (Blue Histogram) or an anti-human MSLN monoclonal antibody (Catalog HC496407, Yellow Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a Goat Anti-Mouse IgG H&L Polyclonal Antibody, FITC (abInScience: MF690414) and cells analysed on a NovoCyt Flow Cytometer.

Detects recombinant MSLN (Catalog No: HC496021) in indirect ELISA.

**Anti-Human MSLN/Mesothelin Antibody (SS1)
binds with MSLN**



Bioactivity