


Product Information	
Catalog No.	FC10998
Product name	Anti-Human ROBO4 Antibody (MAb1)
Applications	ELISA, FCM, Neutralization
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	MAb1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Magic roundabout, ROBO4, Roundabout homolog 4
Accession	Q8WZ75
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



Western blot

Western blot analysis was performed using anti-ROBO4 monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant human ROBO4 (Catalog No: HA259012)
Lane 2: negative control

**Anti-Human ROBO4 Antibody (MAb1)
binds with ROBO4**

Detects recombinant human ROBO4 (Catalog No: HA259012) in indirect
ELISA

