



Product Information	
Catalog No.	FC11011
Product name	Anti-Human PODXL Antibody (Mab47)
Applications	ELISA, FCM, IHC, Neutralization
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	Mab47
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	2.37 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PCLP1, Podocalyxin-like protein 1, GCTM-2 antigen, Gp200, PODXL, Podocalyxin, PC, PCLP, PCLP-1
Accession	O00592
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data

Anti-Human PODXL Antibody (Mab47) binds with PODXL



Detects recombinant human PODXL (Cat No.: HT033012) in indirect ELISA.




Western blot analysis was performed using anti-PODXL monoclonal antibody at 1ug/mL on various samples.

Lane 1: recombinant human PODXL (Catalog No: HT033012)

Lane 2: negative control

Western blot



Flow-cytometry using anti-human PODXL antibody. Hela cells were stained with an irrelevant antibody (Blue Histogram) or an anti-human PODXL monoclonal antibody (Catalog HT033107, Yellow Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, bound antibody was detected using a Goat Anti-Mouse IgG H&L Polyclonal Antibody, FITC (abInScience: MF690414) and cells analysed on a NovoCyte Flow Cytometer.

Flow cytometry