


Product Information	
Catalog No.	FC11015
Product name	Anti-Human CD175/Tn Antibody (SAA1994)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1994
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.92 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Tn antigen saccharide, N-acetyl- α -D-galactosamine, GalNAc α 1-O-Ser/Thr/Tyr, CD175, Tn, Tn-(GalNAc α 1-O-Ser/Thr)
Accession	ChEBI: 40356
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SEC-HPLC detection for Anti-Human CD175/Tn Antibody (ABS1914).

Bioactivity