


Product Information	
Catalog No.	FC11063
Product name	Anti-Mycobacterium tuberculosis pstS1 Antibody (SAA2170)
Applications	ELISA, FCM, SPR
Species reactivity	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SAA2170
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Phosphate-binding protein PstS 1, PBP 1, PhoS1, PstS-1, 38-kDa glycolipoprotein, 38-kDa lipoprotein, P38, Antigen Ag78, Protein antigen B, Pab, pstS1, phoS1, Rv0934, MTCY08D9.05c
Accession	P9WGU1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Mycobacterium tuberculosis pstS1 Antibody (SAA2170)

SDS-PAGE