


Product Information	
Catalog No.	FC11075
Product name	Anti-Human IL11RA Antibody (SAA2202)
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA2202
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL-11RA, sIL-11R, Interleukin-11 receptor subunit alpha, sIL-11RA, IL11RA, IL-11R-alpha, IL-11R subunit alpha, IL-11 receptor subunit alpha, sIL11RA
Accession	Q14626
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human IL11RA Antibody (SAA2202)

SDS-PAGE