

Product Information	
Catalog No.	FC11086
Product name	Anti-Human CLDN4/Claudin-4 Antibody (SAA0509), FITC
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0509
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CPER, CLDN4, Williams-Beuren syndrome chromosomal region 8 protein, CPE-R, CPE-receptor, WBSCR8, Claudin-4, Clostridium perfringens enterotoxin receptor, CPETR1
Accession	O14493
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.