

Product Information	
Catalog No.	FC11096
Product name	Anti-Mouse CD304/NRP1/NRP2 Antibody (YW68.11.26), FITC
Applications	ELISA, FCM
Species reactivity	Human, Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	YW68.11.26
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	NRP1, Neuropilin-1, VEGF165R, NRP, CD304, Vascular endothelial cell growth factor 165 receptor, Vascular endothelial cell growth factor 165 receptor 2, Neuropilin-2, VEGF165R2, NRP2
Accession	P97333, O35375
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.