

Product Information	
Catalog No.	FC11099
Product name	Anti-Human CD337/NCR3 Antibody (BGA1831), FITC
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	BGA1831
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	NKp30, LY117, NCR3, 1C7, NK-p30, CD337, Natural cytotoxicity triggering receptor 3, Activating natural killer receptor p30, Natural killer cell p30-related protein
Accession	O14931
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.